

Accessories for the DC-**STIMULATOR** range

For our DC-**STIMULATOR** devices we offer a wide range of accessories, such as electrodes of various dimensions and shapes, cables and pastes. You can order these products individually or in a set. Please contact us. We will be happy to advise you which electrodes are best suited for your project.



Please send your order by e-mail to: sales@neurocaregroup.com

Content


Accessories for the DC- STIMULATOR range.....	1
Equipment starter set for DC- STIMULATOR / DC- STIMULATOR MOBILE	2
Equipment starter sets for DC- STIMULATOR MC	3
Rubber electrodes standard sizes (rectangular).....	4
Sponge pads for rubber electrodes	5
Value packs of electrodes and sponges	6
Rubber electrodes suitable for MR	7
Rubber electrodes - special dimensions and shapes	8
EEG electrode caps	10
Caps for tDCS and tACS	11
Other materials (rubber straps, cables, pastes ...)	13
Contact for enquiries and orders	14

Equipment starter set for DC-**STIMULATOR** / DC-**STIMULATOR** MOBILE

The equipment starter sets for the DC-**STIMULATOR** and the DC-**STIMULATOR** MOBILE consist of:

- 1 pair of rubber electrodes in different standard sizes (rectangular)
- 1 pair of sponge pads
- 1 pair of connection cables
- 1 combination of elastic rubber straps to attach the electrode to the head

All components can be ordered individually or in value packs (please see the following categories).




equipment starter set: 9 cm²
order # 305018



equipment starter set: 25 cm²
order # 305010



equipment starter set: 35 cm²
order # 305012



equipment starter set: 50 cm²
order # 305014



equipment starter set: 100 cm²
order # 305016

Equipment starter sets for DC-**STIMULATOR** MC

Each equipment starter set for the DC-**STIMULATOR** MC4 consists of:

- 4 pairs of rubber electrodes in different standard sizes (rectangular)
- 4 pairs of sponge pads
- 4 pairs of connection cables
- 4 combinations of elastic rubber straps to attach the electrode to the head

When using the DC-**STIMULATOR** MC8 or MC16 you require 2 or 4 starter sets.

All components can be ordered individually or as value pack (see following categories).



pack of electrodes MC4
3 x 3 cm (9 cm²)
order # 305026




pack of electrodes MC4
5 x 5 cm (25 cm²)
order # 305020



pack of electrodes MC4
5 x 7 cm (35 cm²)
order # 305021

Rubber electrodes standard sizes (rectangular)



3 x 3 cm, 9 cm²
1 pair
305054



5 x 5 cm, 25 cm²
1 pair
305050



5 x 7 cm, 35 cm²
1 pair
305051



5 x 10 cm, 50cm²
1 pair
305052



10x10cm, 100cm²
1 pair
305053

Sponge pads for rubber electrodes

The indicated dimensions are approximate values.



3 x 3 cm, 9 cm², 1 pair
305064



5 x 5 cm, 25cm², 1 pair
305060



5 x 7 cm, 35cm², 1 pair
305061



5 x 10 cm, 50 cm², 1 pair
305062



10x10 cm, 100 cm², 1 pair
305063



3 x 3 cm, 9 cm², 10 pairs
305064-11



5 x 5 cm, 25 cm², 10 pairs
305065



5 x 7 cm, 35 cm², 10 pairs
305066




5 x 10 cm, 50 cm², 10 pairs
305067




10x10 cm, 100 cm², 10 pairs
305068

Value packs of electrodes and sponges




value pack 9 cm²

order # 305112 (pack of 5)
order # 305122 (pack of 10)



value pack 25 cm²

order # 305113 (pack of 5)
order # 305123 (pack of 10)



value pack 35 cm²


order # 305114 (pack of 5)
order # 305124 (pack of 10)

Rubber electrodes suitable for MR

The indicated dimensions are approximate values.



suitable for MR
3 x 3 cm, 9 cm²
1 pair
order # 305059



suitable for MR
5 x 5 cm, 25 cm²
1 pair
order # 305055



suitable for MR
5 x 7 cm, 35 cm²
1 pair
order # 305056



suitable for MR
5 x 10 cm, 50 cm²
1 pair
order # 305057

Rubber electrodes - special dimensions and shapes

The indicated dimensions are approximate values.



Ø: 10 mm
area: 0,8cm²
without hole
2 mm thick
1 pair

order # 305090-01



Ø: 20 mm
area: 3 cm²
without hole
1 mm thick
1 pair

order # 305090-02



Ø: 25 mm
area: 5 cm²
without hole
2 mm thick
1 pair

order # 305090-11



Ø: 34 mm
area: 9 cm²
without hole
2mm thick
1 pair

order # 305090-03



Ø: 45 mm
area: 16 cm²
without hole
2 mm thick
1 pair

order # 305090-04



Ø: 75 mm
area: 44 cm²
without hole
2 mm thick
1 pair

order # 305090-05



Ø: 80 mm
area: 50 cm²
without hole
2 mm thick
1 pair

order # 305090-06



50 mm x 30 mm
area: 37 cm²
1 mm thick
1 pair

order # 305090-14



140mm x 70mm
area: 98 cm²
size: 2 mm
1 pair

order # 305090-10



outer Ø: 45 mm
inner Ø: 15 mm
area: 14 cm²
with hole
2mm thick
1 pair

order # 305090-17



outer Ø: 45 mm
inner Ø: 15 mm
area: 14 cm²
with hole
2mm thick
1 pair

order # 305090-17




outer Ø: 48 mm
inner Ø: 24 mm
area: 15 cm²
with hole
2 mm thick
1 pair

order # 305090-12



outer Ø: 75 mm
inner Ø: 20 mm
area: 40 cm²
with hole
2 mm thick
1 pair

order # 305090-07



outer Ø: 75 mm
inner Ø: 30 mm
area: 37 cm²
with hole
2 mm thick
1 pair

order # 305090-15



outer Ø: 100mm
inner Ø: 70mm
area: 40 cm²
with hole
2 mm thick
1 pair

order # 305090-16



outer Ø: 100mm
inner Ø: 75 mm
area: 40 cm²
with hole
2mm thick
1 pair

order # 305090-08



outer Ø: 110mm
inner Ø: 90mm
area: 31 cm²
with hole
2 mm thick
1pair

order # 305090-13

EEG electrode caps



EEG electrode cap
no slots for adapters
chin strap

order # 905440



EEG electrode cap for ASIA
no slots for adapters
chin strap

order # 905441

Caps for tDCS and tACS

Knitted caps for transcranial electrical stimulation (tES) with direct current (tDCS) or alternating current (tACS) are available in various sizes and configurations.

On the next page you will find all available configurations.



Available sizes:


circumference/ hight*	2/2	2/4	3/3	4/2	4/4
circumference in cm	< 56,5	< 56,5	56,5 - 59,5	> 59,6	> 59,6
hight in cm	> 11,0	> 14,0	11,1 - 13,0	< 11,0	> 14,0

The height is measured with a straight line between the height of the corner of the eye and the highest point on the head.


*Please specify the size when ordering (e.g. "size 3/3").




Available configurations:




Knitted tES cap with chin strap and cables
Electrode positions: C3 and FP2
Size of electrodes: 5 cm x 7 cm
Order number: 905051




Knitted tES cap with chin strap and cables
Electrode positions: C4 and FP1
Size of electrodes: 5 cm x 7 cm
Order number: 905052



Knitted tES cap with chin strap and cables
Electrode positions: F3 and FP2
Size of electrodes: 5 cm x 7 cm
Order number: 905053



Knitted tES cap with chin strap and cables
Electrode positions: F3 and F4
Size of electrodes: 5 cm x 7 cm
Order number: 905054



Knitted tES cap with chin strap and cables
Electrode positions: F3 and F4
Size of electrodes: 5 cm x 5 cm
Order number: 905055

Other materials (rubber straps, cables, pastes ...)



connection cable for rubber electrodes
approx. 150 cm, 1 pair (1-channel)
order # 305069



elastic rubber strap for rubber electrodes, 1 set
order # 305070



Extension rubber strap for rubber electrodes, 1 set
order # 305070-12



Velcro® strip for rubber electrodes, 1 set
order # 305072



extension of connection cables for rubber electrodes (5 m), 1 pair (1-channel)
order # 305075



extension of connection cables for rubber electrodes (10 m), 1 pair (1-channel)
order # 305076



Ten20 electrode paste
3 x 228 g cans
order # 905821



Ten20 electrode paste
3 x 114 g cans
order # 905822

Contact for enquiries and orders

neurocare group AG
Albert-Einstein-Straße 3
98693 Ilmenau

sales@neurocaregroup.com

+49 - 3677 – 68 97 90