

**Checklist for DC-STIMULATOR MC**

	shipping case (borrowable)
--	----------------------------

**PANEL-PC**

	PANEL-PC
	power supply unit
	power cable 3-pole
	keyboard (if available: USB adapter)
	mouse (if available: USB adapter)

**STIMULATOR**

	STIMULATOR/-s
	USB cable to connect STIMULATOR and PANEL-PC
	power cable (red-yellow) between STIMULATOR and PANEL-PC

**Accessories**

	adaptor /-s (black)
	adaptor box/-es
	box cable

**DC-STIMULATOR MC 16 channel**

	cable set to connect two STIMULATOR (TNC cable, BNC cable, power cable)
--	-------------------------------------------------------------------------

**Consumables**

	connecting cable red-blue
	other:

**Accessories DC-STIMULATOR MC MR**

	INNER BOX /-es
	OUTER BOX /-es
	ELECTRODE CABLE
	STIMULATOR CABLE
	BOX CABLE

## Shipping case (borrowable)

Ensuring safe and reliable shipment we provide a borrowable shipping case for the DC-**STIMULATOR** MC (4 and 8 channel). For this please contact the service department.

A manual how to place the components in the shipping case is listed at the end of this document.

## PANEL-PC



## STIMULATOR



## **Power supply**

The power supply consists of the power supply unit and the power cable (3-pole). Please always send both!

power supply unit



power cable (3-pole)



## **Keyboard and mouse**

If available please always enclose the USB-PS2 adaptor

keyboard



mouse



### Connecting cable

power cable (red-yellow)



USB cable




if 16 channel: cable set to connect two STIMULATOR



### Accessories

adaptor box



adaptor



box cable



connecting cable red-blue



Accessories DC-STIMULATOR MC MR

**OUTER BOX**



**INNER BOX**



**STIMULATOR CABLE**



**ELECTRODE CABLE**



**BOX CABLE**



**Packing the DC-STIMULATOR MC in the shipping case**



- 1 STIMULATOR, adaptor/-s plugged in
- 2 adaptor box /-es
- 3 box cable
- 4 power cable (red-yellow)
- 5 USB cable
- 6 power supply unit
- 7 power cable (3-pole)
- 8 mouse
- 9 keyboard



- 10 manual / documents
- 11 PANEL-PC

Please place the connecting cable red-blue and further accessories next to the power supply unit and the mouse.